


# Sample storage



## Tubes, cryogenic vials, 2D barcoded and uncoded, screw cap, sterile, Thermo Scientific Nunc® Cryobank



For cell and protein storage.

- Certified RNase/DNase free, non-pyrogenic
- Automation friendly design
- Unique low cell and protein bind surface on Cryobank vials
- Dense storage format for optimised freezer utilisation
- Unique 14 x 14 array 2D code
- Available 2D coded or uncoded
- CE marked

Unique 14 x 14 array 2D code on tubes, e.g. NU XX 11111 C.

2D code is laser etched into polypropylene bottom of the Cryobank tubes and cannot be lost during storage.

1D code on rack. 96 individual Cryobank tubes with screw caps in a 96 well format frame (SBS footprint). Working volume of 0.5mL or 1.0mL using the same rack. Tubes click into rack allowing one hand operation.

For accessories and spares please contact Customer Service (see contact details on inside of front cover).

Catalogue No	Legacy No	Alt. No	Description	Volume, working, mL	Inner pack qty	Pack qty
10044632	-	374083 <sup>†</sup>	Racked uncoded	0.5	96	192
10126412	-	374084 <sup>†</sup>	Racked 2D coded	1.0	96	192
11702265	-	374086	Racked 2D coded	0.5	96	960
10054632	-	374087	Racked uncoded	0.5	96	960
10136412	-	374088	Racked 2D coded	1.0	96	960
10309822	-	374089	Racked uncoded	1.0	96	960
10257752	-	374080	Uncoded loose pack	0.5	96	960
10320013	-	374081	Uncoded loose pack	1.0	96	960

<sup>†</sup> Includes black manual screwdriver for unscrewing and fastening caps



## Tubes, cryogenic vials, 2D barcoded and uncoded, screw cap, Thermo Scientific Nunc® Bank-It™



- Certified RNase/DNase free, non-pyrogenic
- Automation friendly design
- Low DNA absorption on Bank-It™ vials
- Dense storage format for optimised freezer utilisation
- Unique 14 x 14 array 2D code
- Available with 2D code or uncoded, supplied racked
- CE marked

Unique 14 x 14 array 2D code on tubes, e.g. NU XX 11111 C.

2D code is laser etched into polypropylene bottom of the Bank-It™ tubes and cannot be lost during storage.

1D code on rack. 96 individual Bank-It™ tubes with screw caps in a 96 well format frame (SBS footprint). Working volume of 0.5mL or 1.0mL using the same rack. Tubes click into rack allowing one hand operation.

For accessories and spares please contact Customer Service (see contact details on inside of front cover).

Catalogue No	Legacy No	Alt. No	Description	Volume, working, mL	Inner pack qty	Pack qty
10533933	-	374074	2D coded	0.5	96	960
10298372	-	374075	Uncoded	0.5	96	960
10543933	-	374078	2D coded	1.0	96	960

### Accessories

Catalogue No	Legacy No	Alt. No	Description
10753754	-	374009	Tube selection tool, Bank-It™

