


Spectrophotometers



Continuation

Spectrophotometers, software and accessories, calibration and validation



Calibration

Catalogue No	Legacy No	Alt. No	Description
10224818	-	840-211300	Mercury lamp for validation and calibration, primary standard recommended to test wavelength accuracy of Evolution™ 200 series spectrophotometers
10031518	-	840-213000	Calibration validation carousel with smart base, European Pharmacopeia (EP) compliant, with set of certified reference materials for Evolution™ 200 series spectrophotometers
10771789	-	840-212900	Calibration validation carousel with smart base, US Pharmacopeia (USP) compliant, with set of certified reference materials for Evolution™ 200 series spectrophotometers
11556959	-	10010721	Calibration validation carousel with smart base, NIST traceable, with set of certified reference materials for Evolution™ 200 series spectrophotometers
10040623	-	10010621	Calibration validation carousel with smart base, NPL traceable, with set of certified reference materials for Evolution™ 200 series spectrophotometers

Standards

Catalogue No	Legacy No	Alt. No	Description
12197820	-	9423 185 03001	ND Absorbance filters, set of 4, NIST/NPL traceable
11781161	-	9423 179 38091	ND Absorbance filters, set of 9, NIST/NPL traceable
11791161	-	9423 185 03112	Wavelength filters, holmium and didymium, NIST/NPL traceable
11823654	-	222-226400	UV-Vis wavelength standard, holmium oxide in perchloric acid, NIST traceable

Recalibration and certification


Catalogue No	Legacy No	Alt. No	Description
11853764	-	9423 UV9 1460E	Recalibration service and recertification of calibration validation carousel, for use with CVCs 10010721, 10010701, 9423 UV6 1460E
10213001	-	9423 UV9 1400E	Recalibration service and recertification of calibration validation carousel, for use with CVCs 10010621, 10010601, 9423 UV6 1400E
12993776	-	4013 229 89851	Recalibration service and recertification of 4 neutral density absorbance filters for use with 942318503001
12903786	-	4013 229 89811	Recalibration service and recertification of 4 neutral density absorbance filters for use with 942317938091
11601141	-	4013 229 89841	Recalibration service and recertification of holmium oxide/didymium oxide filters 942318503112

Have you thought about.....

Cuvettes, disposable, PS and PMMA, Fisherbrand®



Catalogue No	Legacy No	Description	Capacity, mL	Material	Inner pack qty	Pack qty
11682599	-	Standard/macro	4.5	Polystyrene	100	500
11602609	-	Semi-micro	1.5	Polystyrene	100	500
11954395	-	4 clear sides	4.5	Polystyrene	100	500
11944385	-	Standard/macro	4.5	Methacrylate	100	500
11904385	-	Semi-micro	1.5	Methacrylate	100	500
11924405	-	4 clear sides	4.5	Methacrylate	100	500



Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.