

Stirrers, magnetic, unheated, C-Mag MS



- Ceramic glass plate, offering excellent chemical resistance
- Powerful motor for quantities up to 15L (H₂O)
- Elevated panel for protection against leaking liquids

Technical Specification - Specific

	10706605	10555824	10002954
Model No	MS 4	MS 7	MS 10
Stirring quantity H ₂ O, L	5	10	15
Speed range, rpm	100 to 1,500	100 to 1,500	100 to 1,500
Plate dimensions, mm	100 x 100	180 x 180	260 x 260
Dimensions [w x d x h], mm	150 x 260 x 105	220 x 330 x 105	300 x 415 x 105
Mass, kg	3	5	6
Electrical supply	230V, 50/60Hz	230V, 50/60Hz	230V, 50/60Hz

Catalogue No	Legacy No	Alt. No	Description
10706605	-	4240200	C-MAG MS 4 stirrer
10555824	-	3487000	C-MAG MS 7 stirrer
10002954	-	4240400	C-MAG MS 10 stirrer

Accessory

Catalogue No	Legacy No	Alt. No	Description	Dimensions, mm
10321352	-	1293100	Stir bar remover	350 x 8 (l x dia.)

Have you thought about.....

Magnetic followers, PTFE, value range, Fisherbrand®



Catalogue No	Legacy No	Length, mm	Diameter, mm	Pack qty
11587572	-	12	3	10
11597572	-	25	8	10
11507582	-	40	8	10
11517582	-	50	8	10



Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.