

Magnetic followers, PTFE, value range, Fisherbrand®


Fisherbrand®



| Catalogue No | Legacy No | Length, mm | Diameter, mm | Pack qty |
|--------------|-----------|------------|--------------|----------|
| 11587572 | - | 12 | 3 | 10 |
| 11597572 | - | 25 | 8 | 10 |
| 11507582 | - | 40 | 8 | 10 |
| 11517582 | - | 50 | 8 | 10 |

Magnetic followers, glass, Fisherbrand®

Fisherbrand®



| Catalogue No | Legacy No | Length, mm | Diameter, mm | Pack qty |
|--------------|-----------|------------|--------------|----------|
| 11537792 | - | 6 | 5 | 5 |
| 11547792 | - | 12 | 5 | 5 |
| 11557792 | - | 19 | 5 | 5 |
| 11567792 | - | 25 | 6 | 5 |
| 11587792 | - | 50 | 8 | 5 |

Magnetic followers, PTFE, centre rim, Fisherbrand®

Fisherbrand®



| Catalogue No | Legacy No | Length, mm | Diameter, mm | Pack qty |
|--------------|-----------|------------|--------------|----------|
| 11587802 | - | 12 | 4.5 | 5 |
| 11597802 | - | 20 | 6 | 5 |
| 11507812 | - | 25 | 6 | 5 |
| 11517812 | - | 35 | 6 | 5 |
| 11527812 | - | 40 | 8 | 5 |
| 11537812 | - | 45 | 8 | 5 |
| 12641156 | - | 50 | 8 | 5 |
| 10018131 | - | 60 | 10 | 1 |
| 10373592 | - | 70 | 10 | 1 |

Magnetic followers, PTFE, octagonal, Fisherbrand®

Fisherbrand®



The flat surfaces provide a greater mixing action. All have pivot rings. AlNiCo-V magnets encased in PTFE.

| Catalogue No | Legacy No | Length, mm | Diameter, mm | Pack qty |
|--------------|-----------|------------|--------------|----------|
| 11858872 | - | 13 | 3 | 10 |
| 11577772 | - | 13 | 8 | 10 |
| 11923574 | - | 15 | 8 | 10 |
| 11943574 | - | 25 | 8 | 10 |
| 11507782 | - | 38 | 8 | 10 |
| 11517782 | - | 51 | 8 | 10 |
| 11527782 | - | 25 | 9.5 | 10 |
| 11537782 | - | 38 | 9.5 | 10 |
| 11547782 | - | 51 | 9.5 | 10 |
| 11557782 | - | 64 | 9.5 | 10 |

Magnetic followers, PTFE, triangular, Fisherbrand®

Fisherbrand®




| Catalogue No | Legacy No | Length, mm | Diameter, mm | Pack qty |
|--------------|-----------|------------|--------------|----------|
| 11577812 | - | 12 | 6 | 5 |
| 11587812 | - | 25 | 8 | 5 |
| 11597812 | - | 35 | 9 | 5 |
| 11848872 | - | 50 | 12 | 5 |

Magnetic followers, PTFE, wing, Fisherbrand®

A

Fisherbrand®

For effective stirring in small tapered and round bottom test, centrifuge and sample tubes.



| Catalogue No | Legacy No | Height, mm | Diameter, mm | Tube diameter [maximum], mm |
|--------------|-----------|------------|--------------|-----------------------------|
| 10042853 | - | 9 | 5.5 | 6 to 7 |
| 10236443 | - | 11 | 8 | 9 to 10 |
| 10356262 | - | 16 | 10 | 11 to 12 |
| 10653283 | - | 10 | 13 | 14 to 15 |
| 10526489 | - | 20 | 13 | 14 to 15 |