

Magnetic followers, PTFE, single sided star head, Fisherbrand®

Fisherbrand®




Catalogue No	Legacy No	Length, mm	Diameter, mm	Pack qty
11527822	-	8	10	5
11537822	-	12	14	5
11547822	-	13	17	5
11557822	-	15	25	5
10690443	-	17	40	1
10600453	-	17	60	1

Magnetic followers, PTFE, micro, Fisherbrand®

Fisherbrand®

For stirring in test tubes, small size beakers or other containers.

Supplied in re-useable snap cap vials.




Catalogue No	Legacy No	Length, mm	Diameter, mm	Pack qty
11858882	-	2	2	10
11878882	-	5	2	10
11898882	-	8	1.5	10
11888952	-	7	2	10
11318703	-	10	3	10

Magnetic followers, PTFE, Saturn® Spinbar, Fisherbrand®



Fisherbrand®

Manufactured from inert and chemically resistant PTFE the Saturn® Spinbar is autoclavable at 121°C. The Saturn® has a large bulb at the centre which lifts the stir bar at the ends whilst spinning, and thereby reduces the surface contact area allowing the bar to spin with decreased drag. Available in 2 lengths to match a wide variety of vessels, with both round and flat bottoms.



Catalogue No	Legacy No	Length, mm
11567572	-	40
11577572	-	50

Magnetic followers, PTFE, cross, Fisherbrand®



Fisherbrand®

A stable, general purpose stirrer. Provides a deep, quiet vortex even at relatively high speeds.



Catalogue No	Legacy No	Length, mm	Diameter, mm	Pack qty
11521983	-	10	5	5
11808942	-	20	8	5
11818942	-	25	9	5
11828942	-	30	10	5
11838942	-	38	11	5

Magnetic followers, PTFE, oval, Fisherbrand®

Fisherbrand®



Catalogue No	Legacy No	Length, mm	Diameter, mm	Pack qty
11597792	-	10	5	5
11507802	-	15	6	5
11517802	-	20	10	5
11527802	-	25	12	5
11338703	-	30	16	5
11818882	-	35	16	5
10603731	-	40	20	1
10664481	-	50	20	1
10055981	-	64	20	1
10609162	-	70	20	1

Magnetic followers, PTFE, coloured, double ended, Fisherbrand®



Fisherbrand®

Ideal for colour coded experiments. The double paddle action gives efficient stirring plus high stability.



Note: Coloured PTFE coatings are not as inert as pure (white) PTFE

Catalogue No	Legacy No	Length, mm	Diameter, mm	Colour
10799681	-	35	8	White
10877880	-	35	8	Red
10736690	-	35	8	Blue
10062719	-	35	8	Yellow
10271861	-	55	8	White
10586319	-	55	8	Red
10144189	-	55	8	Blue
10746690	-	55	8	Yellow