

Magnetic followers, PTFE, octagonal, colour coded, Fisherbrand®

Fisherbrand®



Colour coded, octagonal stirrer bars. For use where identification is of prime importance. Note: Coloured PTFE coatings are not as inert as pure PTFE.

Standard sizes

Catalogue No	Legacy No	Description	Pack qty
11888902	-	Red, 13mm x 8mm	10
11898902	-	Red, 15mm x 8mm	10
11572303	-	Red, 25mm x 8mm	10
11582303	-	Red, 38mm x 8mm	10
11792145	-	Blue, 13mm x 8mm	10
12574987	-	Blue, 15mm x 8mm	10
11512313	-	Blue, 25mm x 8mm	10
11522313	-	Blue, 38mm x 8mm	10
11532313	-	Blue, 51mm x 8mm	10
11542313	-	Yellow, 13mm x 8mm	10
11898912	-	Yellow, 15mm x 8mm	10
11562313	-	Yellow, 25mm x 8mm	10
11572313	-	Yellow, 38mm x 8mm	10
11582313	-	Yellow, 51mm x 8mm	10

Micro sizes

Catalogue No	Legacy No	Description	Pack qty
11507183	-	Red, 2mm x 2mm	10
12696316	-	Red, 5mm x 2mm	10
11527183	-	Red, 7mm x 2mm	10
11537183	-	Red, 8mm x 1.5mm	10
11547183	-	Red, 10mm x 3mm	10
11557183	-	Red, 15mm x 1.5mm	10
11567183	-	Blue, 2mm x 2mm	10
11577183	-	Blue, 5mm x 2mm	10
11587183	-	Blue, 7mm x 2mm	10
11597183	-	Blue, 8mm x 1.5mm	10
12686316	-	Blue, 10mm x 3mm	10
11517193	-	Blue, 15mm x 1.5mm	10
11527193	-	Yellow, 2mm x 2mm	10
11537193	-	Yellow, 5mm x 2mm	10
11547193	-	Yellow, 7mm x 2mm	10
11557193	-	Yellow, 10mm x 3mm	10
11567193	-	Yellow, 15mm x 1.5mm	10

Magnetic followers, rare earth samarium cobalt, Fisherbrand®



Fisherbrand®



Stir bars using rare earth magnets have significantly increased magnetic strength compared to those using AlNiCo magnets.

They are almost totally resistant to demagnetisation. Stir bars using rare earth magnets are identified by a chemically inert, carbon black spot. N.B. The high magnetic coupling power of rare earth stir bars may cause abrasion of the PTFE coating against the base of the container and where stirring arrays are used, as in multi-well plates, lateral interaction between stir bars may be observed.

Elliptical

Catalogue No	Legacy No	Length, mm	Diameter, mm	Pack qty
11501953	-	10	6	5
11888922	-	15	10	5
11848922	-	25	14	5
10341542	-	50	24	1

Cylindrical

Catalogue No	Legacy No	Length, mm	Diameter, mm	Pack qty
11541953	-	8	3	5
11551953	-	12	4.5	5
11561953	-	25	6	5
10582972	-	40	8	5
10456402	-	60	10	1

Magnetic followers, set, PTFE, octagonal, Fisherbrand®

Fisherbrand®

Set of 12 assorted octagonal followers with pivot ring, in a compartmented box. Set comprises two each of the following sizes:
113mm x 8mm
16mm x 8mm
25mm x 10mm
38mm x 10mm
51mm x 10mm
64mm x 10mm



Catalogue No	Legacy No	Description
11547772	-	Set of 12 magnetic followers, octagonal