

Test tubes, borosilicate glass, round bottom, screw thread with polypropylene caps, disposable, Fisherbrand®



Ideal for tissue culture, bacteriology, clinical chemistry, blood typing and cross matching procedures.

With natural, linerless, polypropylene screw cap and marking spot. Tubes packed in shrink wrapped trays with caps packed separately.

Catalogue No	Legacy No	Diameter [external], mm	Length, mm	GPI thread finish	Pack qty
11517413	-	13	100	13-415	1,000
11527413	-	16	100	15-415	1,000
10421541	-	16	125	15-415	1,000
11557413	-	16	150	15-415	1,000
11567413	-	20	125	18-415	500
11577413	-	20	150	18-415	500



Test tubes, borosilicate glass, round bottom, screw thread without caps, disposable, Fisherbrand®



Ideal for tissue culture, bacteriology, clinical chemistry, blood typing and cross matching procedures.

Available with or without marking spot. Tubes packed in shrink wrapped trays.

With marking spot

Catalogue No	Legacy No	Diameter [external], mm	Length, mm	GPI thread finish	Pack qty
11587413	-	13	100	13-415	1,000
12327279	-	16	100	15-415	1,000
11527423	-	16	125	15-415	1,000
11537423	-	16	150	15-415	1,000
11547423	-	20	125	18-415	500
11557423	-	20	150	18-415	500

Without marking spot

Catalogue No	Legacy No	Diameter [external], mm	Length, mm	GPI thread finish	Pack qty
11567423	-	13	100	13-415	1,000
11577423	-	16	150	15-415	1,000
11587423	-	20	150	18-415	500



Test tubes, borosilicate glass, light walled, rimless, Fisherbrand®



Premium quality tubes with sturdy, uniform bottoms and consistent lengths. The 10mm x 75mm and 12mm x 75mm sizes are suitable for cell washing procedures. The 6mm x 50mm size is often referred to as a 'Durham tube'. Approximate wall thickness is 0.6mm. 10022253 is manufactured from flint glass.

Catalogue No	Legacy No	Diameter [external], mm	Length, mm	Pack qty
10022253	-	6	50	1,000
12347279	-	10	75	1,000
11517403	-	12	75	1,000
11527403	-	13	100	1,000
11537403	-	15	85	1,000
11547403	-	16	100	1,000
11557403	-	16	125	1,000
11577403	-	16	150	1,000
11587403	-	18	150	500
11597403	-	20	150	500
11507413	-	25	150	500



For Fisherbrand® closures and caps, please see pages 1541 to 1542

