



Vials, sample, clear and amber, Type 1 Class A borosilicate glass, with fitted caps, Fisherbrand®



- Supplied in partitioned trays
- Standard vials have fitted black phenolic caps with either polycone liner or PTFE faced white rubber liner
- Also available: pre-cleaned EPA (Environmental Protection Agency) vials for VOA (volatile organic analysis), with or without certification, having fitted white polypropylene open top caps with PTFE faced silicone septa

Clear sample vials with fitted black phenolic polycone lined caps

Catalogue No	Legacy No	Capacity, mL	Capacity, drams	O.D., mm	Height, mm	GPI thread finish	Pack qty
11660112	-	4	1	15	48	13-425	144
11537733	-	8	2	17	63	15-425	144
11547733	-	12	3	19	68	15-425	144
11553522	-	16	4	21	73	18-400	144
11583522	-	20	5	28	60	24-400	72
11563522	-	24	6	23	88	20-400	144
11593522	-	25	6.25	28	73	24-400	72
11573522	-	30	8	25	94	22-400	144

Clear sample vials with fitted black phenolic PTFE faced rubber lined caps

Catalogue No	Legacy No	Capacity, mL	Capacity, drams	O.D., mm	Height, mm	GPI thread finish	Pack qty
11503532	-	2	0.5	12	38	8-425	288
11523532	-	4	1	15	48	13-425	144
11543532	-	8	2	17	63	15-425	144
11563532	-	12	3	19	68	15-425	144
11573532	-	16	4	21	73	18-400	144
11513542	-	20	5	28	60	24-400	72
11583532	-	24	6	23	88	20-400	144
11533542	-	25	6.25	28	73	24-400	72
11593532	-	30	8	25	94	22-400	144

Clear short sample vials with fitted black phenolic PTFE faced rubber lined caps

Catalogue No	Legacy No	Capacity, mL	Capacity, drams	O.D., mm	Height, mm	GPI thread finish	Pack qty
11543542	-	2	0.5	15	31	13-425	200
11553542	-	4	1	17	41	15-425	200
11563542	-	6	1.5	19	43	15-425	200

Amber sample vials with fitted black phenolic PTFE faced rubber lined caps

Catalogue No	Legacy No	Capacity, mL	Capacity, drams	O.D., mm	Height, mm	GPI thread finish	Pack qty
11573542	-	2	0.5	12	38	8-425	288
11583542	-	4	1	15	48	13-425	144
11309493	-	8	2	17	63	15-425	144
11503552	-	20	5	28	60	24-400	72
11513552	-	40	10	28	98	24-400	72

Pre-cleaned EPA (Environmental Protection Agency) vials for VOA in both clear and amber, with fitted white polypropylene open top cap and fitted PTFE faced silicone septa. All EPA vials have an outside diameter of 28mm, 24-400 GPI thread finish.

Without certification

Catalogue No	Legacy No	Capacity, mL	Capacity, drams	Colour	Height, mm	Pack qty
11553552	-	40	10	Clear	98	144
11563552	-	40	10	Amber	98	144

With certification

Catalogue No	Legacy No	Capacity, mL	Capacity, drams	Colour	Height, mm	Pack qty
11583552	-	40	10	Clear	98	144
11593552	-	40	10	Amber	98	144