

Vials, clear, with rubber lined caps



Solid top black phenolic caps with 14B white rubber liners. Vials are manufactured from Wheaton-33® low extractable borosilicate glass (meets ASTM Type 1 Class A and USP Type 1 standards) to protect against changes in pH and maintain purity of contents.



Catalogue No	Legacy No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
10242431	-	224881	2	12	38	8	425	288
10739151	-	224882	4	15	48	13	425	144
10160861	-	224884	8	17	63	15	425	144
11734159	-	224885	12	19	68	15	425	144
10504291	-	224886	16	21	73	18	400	144

Vials, amber, with rubber lined caps



Solid top black phenolic caps with 14B white rubber liners. Vials are manufactured from Wheaton '320' amber glass (meets ASTM Type 1 Class B and USP Type 1 standards) to protect light-sensitive samples and maintain purity of contents.



Catalogue No	Legacy No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
10122531	-	224981	2	12	38	8	425	288
10514481	-	224982	4	15	48	13	425	144
10612332	-	224984	8	17	63	15	425	144

Vial caps, black phenolic, screw fit



Black phenolic screw caps with PTFE faced rubber liners for optimum chemical resistance. Do not autoclave containers with caps attached tightly.

Catalogue No	Legacy No	Alt. No	Closure size, mm	GPI thread finish	Pack qty
11511235	-	240406	8	425	200
10778571	-	240408	13	425	200
10058860	-	240409	15	425	200
10246862	-	240415	18	400	200
10742973	-	240418	24	400	100



Serving Science through our catalogue ranges:

- Chemicals
- Laboratory
- Life Science
- Safety



ONE SOURCE. INFINITE SOLUTIONS.®