

## Vials, borosilicate glass, liquid scintillation, with caps, 20mL, Fisherbrand®



**Fisherbrand®**

Potassium-free 33 borosilicate glass provides consistently low activity counts and excellent light transmission. Fit all common counters. Caps fit GPI 22-400 thread finish, and are suitable for marking. Vials with caps attached are shrinkwrapped in five trays of 100 per case. When caps are separate, vials are shrinkwrapped in five trays of 100 per case, and caps are packed in five polybags of 100 caps each.

Vial dimensions, with cap (O.D. x h): 28mm x 61mm.

Catalogue No	Legacy No	Cap material	Cap style	Packaging type	Pack qty
12353317	-	Urea	Cork-backed metal foil liner	Caps attached	500
12363317	-	Urea	Cork-backed metal foil liner	Caps separate	500
12373317	-	Urea	Cone-shaped plastic liner	Caps separate	500
12383317	-	Polypropylene	Pulp-backed metal foil liner	Caps attached	500
12393317	-	Polypropylene	Pulp-backed metal foil liner	Caps separate	500
12303327	-	Polyethylene	No liner	Caps separate	500



## Vials, borosilicate glass, scintillation, 20mL, Thermo Scientific Chromacol®

**Thermo  
SCIENTIFIC**

- Noise = 2.28
- Background count = 13CPM
- Quenching Index Factor = 49

Vial with high reflection foil lined caps

Catalogue No	Legacy No	Alt. No	Capacity, mL	Colour	Diameter, mm	Height, mm	Pack qty
11357371	-	20-EPSCA	20	Clear	27	57	500

Spare caps, white, urea

Catalogue No	Legacy No	Alt. No	Description	Pack qty
11543690	-	20-EPSCA	High reflection foil lined	500
12316928	-	20-EPSP	Polyethylene lined	500



## Vials, glass, liquid scintillation

**WHEATON**

The original scintillation vial, 20mL capacity, manufactured from Wheaton '180' low potassium glass under critical manufacturing tolerances. Background counts are consistent and low whilst UV transmission is high. Tight fitting white screw caps supplied. Dimensions 28mm x 61mm (with cap on).



Catalogue No	Legacy No	Alt. No	Closure size, mm	Closure	Cap style	Inner pack qty	Pack qty
11748819	-	986541	22	Polypropylene	Metal foil	100	500
11758819	-	986542	22	Urea	Metal foil	100	500
10458043	-	986546	22	Urea	Poly-seal cone	100	500
10006761	-	986548	22	Urea	Polyethylene disc	100	500
11768819	-	986561	24	Polypropylene	Metal foil	100	500
11506325	-	986568	24	Urea	Polyethylene	100	500