

Water purification systems, Type 2, Thermo Scientific Barnstead®, Pacific™ TII



Accessories

Catalogue No	Legacy No	Alt. No	Description
11327151	-	09.4000	Pre-treatment complete with a 5µm filter with activated carbon cartridge and a hardness stabiliser cartridge
11357151	-	09.4002	Replacement UV lamp for Pacific™ TII, LabTower™ TII and LabTower™ RO
11317161	-	09.5002	Replacement UV lamp for tank UV housing no. 06.5006
11835912	-	9.2202	Disinfection solution, tablet form, 12 pieces
11397151	-	09.4011	Filter set
11347101	-	06.5201	Replacement combi cartridge (5µm filter and activated carbon), 10 inch
11387101	-	06.5452	Hardness stabilising unit, 10 inch (for pre-treatment no. 09.4000)
11336991	-	19.0046	Pressure pump for 3L/hr to 20L/hr units
11346991	-	19.0050	Pressure pump for 40L/hr units
11337161	-	22.0046	RO membrane for units with permeate performance of 20L/hr or less (15L/hr and 20L/hr units require two)
11377161	-	22.0087	RO membrane for units with permeate performance of 21L/hr to 40L/hr (30L/hr and 40L/hr units require two)
11367001	-	06.1250	Mix Multi Mini, capacity at 10° dH 2.000L
11396981	-	06.1000	Hardness measuring kit
11377081	-	06.5033	Storage tank 30L, as bench top, with level display
11337091	-	06.5063	Storage tank 60L, as bench top, with level display
11397091	-	06.5083	Storage tank 100L, as bench top, with level display
11367081	-	06.5032	Storage tank 30L as 06.5033 but with recirculation/pressure pump
11327091	-	06.5062	Storage tank 60L as 06.5063 but with recirculation/pressure pump
11387091	-	06.5082	Storage tank 100L as 06.5083 but with recirculation/pressure pump
11367061	-	06.5001	Sterile overflow for storage tank 30L to 100L
11337071	-	06.5002	CO ₂ absorber and vent filter for storage tank 30L to 100L
11367071	-	06.5003	Sterile vent filter for storage tank 30L to 100L
11317081	-	06.5006	UV light for tank 30L to 100L
11327081	-	06.5015	Wall mounting for storage tank (30L)
11357081	-	06.5016	Wall mounting for storage tank (60L)

IQ/OQ Qualification documentation available, please enquire.


Have you thought about.....

Cylinders, borosilicate glass, spouted, Class A, blue ceramic graduations, Fisherbrand®



Graduated in blue ceramic. Hexagonal base. Individual reorder code on each item. DIN 12680 BS 604 ISO 4788

Catalogue No	Legacy No	Capacity, mL	Graduations, mL	Pack qty
12952310	-	5	0.1	2
11517832	-	10	0.2	2
11527832	-	25	0.5	2
11537832	-	50	1	2
11547832	-	100	1	2
11557832	-	250	2	2
11577832	-	1,000	10	2
12962320	-	2,000	20	2



Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.