

Test strips, water hardness, Aquadur



Catalogue No	Legacy No	Alt. No	Description	Graduations	Pack qty
10632381	-	91201	Water hardness high	<3, >5, >10, >15, >20, >25 d	100
10028220	-	91220	Water hardness low	<3, >4, >7, >14, >21 d	100



Test sticks, Quantofix®



- Rapid: most tests can be carried out in 10s to 120s
- Convenient: ready to use kits, pre-calibrated, contain all necessary equipment and reagents
- Precise: colour charts adjusted and checked using certified standard solutions
- Shelf life: 2.5 years, store below 30°C



Catalogue No	Legacy No	Alt. No	Description	Graduations	Pack type	Pack qty
12856880	-	91307	Aluminium	0, 5, 20, 50, 200, 500	Sticks and reagents	100
12792771	-	91315	Ammonium	0, 10, 25, 50, 100, 200, 400	Sticks and reagents	100
11874492	-	91332	Arsenic 50	0, 0.05, 0.1, 0.5, 1, 1.7, 3	Sticks and reagents	100
10581363	-	91334	Arsenic 10	0, 0.01, 0.025, 0.05, 0.1, 0.5	Sticks and reagents	100
10474491	-	91314	Ascorbic acid	0, 50, 100, 200, 300, 500, 1,000, 2,000	Sticks	100
10685811	-	91324	Calcium	0, 10, 25, 50, 100	Sticks and reagents	60
10675811	-	91323	Carbonate	0, 3, 6, 10, 15, 20d	Sticks	100
10167430	-	91321	Chloride	0, 500, 1,000, 1,500, 2,000, >3,000	Sticks	100
10100441	-	91301	Chromate	0, 3, 10, 30, 100	Sticks and reagents	100
10490611	-	91303	Cobalt	0, 10, 25, 50, 100, 250, 500, 1,000	Sticks	100
10536732	-	91304	Copper	0, 10, 30, 100, 300	Sticks	100
12826820	-	91318	Cyanide	0, 1, 3, 10, 30	Sticks and reagents	100
10554598	-	91335	EDTA	0, 100, 200, 300, 400	Sticks	100
10238390	-	91328	Formaldehyde	0, 10, 20, 40, 60, 100, 200	Sticks and reagents	100
12843961	-	91336	LubriCheck	0, 15, 50, 75, 130, 200mmol/L KOH	Sticks	100
10063171	-	91325	Molybdenum	0, 5, 20, 50, 100, 250	Sticks and reagents	100
10100491	-	91305	Nickel	0, 10, 25, 50, 100, 250, 500, 1,000	Sticks	100
10384071	-	91313	Nitrite/nitrate	0, 10, 25, 50, 100, 250, 500 for NO ₃ ⁻ 0, 1, 5, 10, 20, 40, 80 for NO ₂ ⁻	Sticks	100
10474841	-	91311	Nitrite	0, 1, 5, 10, 20, 40, 80	Sticks	100
10632321	-	91338	Nitrite pH	0, 1, 5, 10, 20, 40, 80 for NO ₂ ⁻ pH 6.0, 6.4, 6.7, 7.0, 7.3, 7.6, 7.9, 8.2, 8.4, 8.6, 8.8, 9.0, 9.3, 9.6	Sticks	100
10127480	-	91322	Nitrite 3000	0, 0.1, 0.3, 0.6, 1, 2, 3g/L	Sticks	100
10233952	-	91333	Peroxide 1000	0, 50, 150, 300, 500, 800, 1,000	Sticks	100
10159960	-	91312	Peroxide 100	0, 1, 3, 10, 30, 100	Sticks	100
10528482	-	91319	Peroxide 25	0, 0.5, 2, 5, 10, 25	Sticks	100
12773884	-	91320	Phosphate	0, 3, 10, 25, 50, 100 PO ₄ ³⁻	Sticks and reagents	100
10696351	-	91316	Potassium	0, 200, 400, 700, 1,000, 1,500	Sticks and reagents	100
12752613	-	91337	Quat	0, 10, 25, 50, 100, 250, 500, 1,000	Sticks	100
10278340	-	91329	Sulfate	<200, >400, >800, >1,200, >1,600	Sticks	100
10271611	-	91306	Sulfite	0, 10, 25, 50, 100, 250, 500, 1,000	Sticks	100
10771801	-	91309	Tin	0, 10, 25, 50, 100, 250, 500	Sticks	100
10411591	-	91310	Zinc	0, 2, 5, 10, 25, 30, 100	Sticks and reagents	100

For pH meters,
electrodes, buffers
and indicator papers,
please see pages 548
to 607 of our new
Electrochemistry section.

