



### Wipes, tub, pre-saturated alcohol, non-sterile

04

Ideal for use in clean areas, laboratories and semi controlled environments.

Low lint hardwearing tub wipes, pre-saturated with 70/30 isopropanol alcohol, 24g viscose fabric. Fast acting product to sanitise and clean sensitive areas and products.

Catalogue No	Description
CMH-610-060M	Tub, 200 wipes



### Wipes, paper, KIMTECH Science™, Delicate Task

02



- Low lint
- 2-ply

Catalogue No	Dimensions, mm	Sheets per box	Pack qty
CMC-650-052K	200 x 213	200	24



### Wipes, cleanroom, Class 100, ISO Class 5

04



Durx® 770

Suitable for critical cleaning and wiping. Particularly suitable for heavy-duty use, regular maintenance tasks, photo resist applications and around etching baths. Non-woven cellulose and polyester wipers, specially processed to reduce particle generation and increase absorbency while retaining high wet strength.

Durx® 670

Suitable for cleaning diffusion furnaces and as a quartz-ware wiper in implantation areas. Non-woven cellulose and polyester, with low particle generation, high absorbency and high wet strength.

Catalogue No	Type	Dimensions, mm	Pack qty
CMH-600-010L	DURX 770	230 x 230	300
CMH-600-020Y	DURX 670	230 x 230	300



### Wipes, cleanroom, Class 10, ISO Class 4, Gamma Wipe® 120

04



**S**

100% continuous filament, knitted non-run polyester.

Laundered in a Class 10 cleanroom prior to sterilisation. Incredibly strong for tough cleaning. Low particle generation for critical wiping.

Catalogue No	Dimensions, mm	Pack qty
CMH-600-075G	300 x 300	20



### Wipes, cleanroom, Class 100, ISO Class 5, Gamma Wipe®

04



**S**

Non-woven blend of polyester and cellulose.

Cleanliness, durability, and absorbency make it ideal for general sterile clean room applications. Extremely absorbent.

Catalogue No	Dimensions, mm	Pack qty
SAT-600-080R	300 x 300	20

