



Contamination Control

Hand Protection

Glove liners, cleanroom

03




100% continuous filament knitted polyester, designed to be worn under latex, PVC or other outer gloves for maximum comfort and dexterity. One size. May be cleanroom laundered.

Catalogue No	Size	Type	Pack type	Pack qty
SAR-640-010H	Standard	Half-fingered	Pair	10
SAR-640-020E	Large	Half-fingered	Pair	10
SAR-640-040V	Standard	Full-fingered	Pair	10
SAR-640-050S	Large	Full-fingered	Pair	10

Gloves, vinyl, powder free, non-sterile, Profile 3000™ EN 455

50



Protein-free examination glove. Rolled cuff to prevent slippage. Ideal for use in clean areas, laboratories and semi controlled environments. Exceptional grip when both wet and dry. Ambidextrous.

Catalogue No	Size	Pack qty
CMK-195-040L	Small	100
CMK-195-050Y	Medium	100
CMK-195-060F	Large	100
CMK-195-070C	XL	100

Gloves, vinyl, powder free, non-sterile, BioClean ISO Class 4, cleanroom compatible EN 455

04



Fully compatible with Class 10 operating environments. Double bagged in anti-static bags for critical environment use, they are washed in de-ionised water and packed under Class 10 conditions. Extremely low levels of particulates. Manufactured from polyvinyl chloride they are latex-free. Formulated for comfort, excellent dexterity and superior strength. Micro-roughened for improved grip. Ambidextrous. Length 300mm.

Catalogue No	Size	Pack qty
CMK-190-050P	Medium	100
CMK-190-060M	Large	100
CMK-190-070J	XL	100

Gloves, cleanroom, vinyl

04



Powder free, ambidextrous, seamless, disposable vinyl gloves. They feature a rolled cuff to facilitate donning and removal. Available in two grades, 'powder-free' and 'Class 100'. The latter are specially treated and packed in a cleanroom 100 facility to reduce total particulate count.

Clear, anti-static, standard length

Catalogue No	Description	Grade	Size	Pack qty
CMK-870-010Y	Standard	Powder free	Small	100
CMK-870-020F	Standard	Powder free	Medium	100
CMK-870-030C	Standard	Powder free	Large	100

Pink, vinyl, long length (305mm) anti-static gloves. Made to Class 100 cleanroom specifications

Catalogue No	Description	Grade	Size	Pack qty
STP-210-010J	Long	Class 100	Small	100
STP-210-020G	Long	Class 100	Medium	100
STP-210-030D	Long	Class 100	Large	100