



## Contamination Control Hand Protection

Gloves, powder free, BioClean, ambidextrous,  
EN 455 Class/Medical

04



Medical grade examination glove giving extra length to ensure positive seal at wrist. Rolled cuff to prevent slipping. Suitable for clean areas, laboratories and semi controlled environments. Manufactured to exceed requirement of European Glove Standard EN 455 (1 & 2). Tested to AQL 1.5 level. Extractable latex protein contents have been tested in accordance with the ASTM Modified Lowry Test Method D 5712-95 and are less than 50µg. Ambidextrous. Length 305mm.

Catalogue No	Size	Pack qty
CMK-105-030Y	Small	100
CMK-105-021J	Medium	100
CMK-105-040F	Large	100

Gloves, polyisoprene, sterile, hand  
specific, EN455 Class 1 Medical

S3



Made from synthetic polyisoprene and hand specific for maximum comfort. Micro roughened surface for optimum grip. Gamma irradiated, AQL 1.0, meets ASTM F1671-97 Viral Penetration Test.

Catalogue No	Size	Pack type	Pack qty
SAR-681-020Q	5.5	Pair	40
SAR-681-030N	6	Pair	40
SAR-681-040K	6.5	Pair	40
SAR-681-050H	7	Pair	40
SAR-681-060E	7.5	Pair	40
SAR-681-070B	8	Pair	40
SAR-681-080V	8.5	Pair	40
SAR-681-090S	9	Pair	40

Gloves, polychloroprene, Sterile,  
Class 100, ISO Class 5, hand specific  
PPE Category III

S3



The ideal choice for handling cytostatics.

The DermaShield® offers a massive reduction of potential skin reactions and allergies.

It minimises the risks of Type IV allergies because it is free of accelerators, and Type I allergies because it is free of natural rubber. Its polychloroprene composition offers a high level of resistance to ethanol and IPA.

Catalogue No	Size	Pack type	Pack qty
SAR-335-020F	6	Pair	200
SAR-335-030C	6.5	Pair	200
SAR-335-040W	7	Pair	200
SAR-335-050T	7.5	Pair	200
SAR-335-060Q	8	Pair	200
SAR-335-070N	8.5	Pair	200
SAR-335-080K	9	Pair	200

Gloves, polychloroprene, hand specific,  
PPE Category III

S3



The DermaShield® offers a massive reduction of potential skin reactions and allergies.

Minimises the risks of Type IV allergies. Free of accelerators and Type I allergies and free of natural rubber. Its polychloroprene composition offers a high level of resistance to ethanol and IPA. It is the ideal choice for handling cytostatics. Beaded cuff, 310mm long.

Catalogue No	Size	Pack type	Pack qty
SAR-340-020F	6	Pair	200
SAR-340-030C	6.5	Pair	200
SAR-340-040W	7	Pair	200
SAR-340-050T	7.5	Pair	200
SAR-340-060Q	8	Pair	200
SAR-340-070N	8.5	Pair	200
SAR-340-080K	9	Pair	200