



Fisherbrand®
Cylinders, borosilicate glass, spouted, hexagonal base, left handed, BS 604, ISO4788

- Graduated in blue ceramic for visibility when pouring with left hand
- Individual re-order code on each product
- Supplied in packs of 2

Cylinders, Pyrex® glass, spouted, detachable base

- White graduations
- Cylinder body and base available separately
- Supplied in packs of 2



Catalogue No	Cap., mL	Grads., mL
FB56441	5	0.1
FB56443	10	0.2
FB56445	25	0.5
FB56447	50	1.0
FB56449	100	1.0
FB56451	250	2.0
FB56453	500	5.0
FB56455	1,000	10.0
FB56457	2,000	20.0

Catalogue No	Cap, mL	Grads., mL
CYP-330-070X	50	1.0
CYP-330-090Y	100	1.0
CYP-330-110F	250	2.0
Glass cylinders only		
CYP-371-070D	50	1.0
CYP-371-090U	100	1.0
CYP-371-110R	250	2.0
Plastic base, pack of 10		
CYP-415-071F	50	-
CYP-415-091C	100	-
CYP-415-111L	250	-

Cylinders, polypropylene, squat form, self draining base

- Squat design for increased stability
- Unique self-draining base eliminating water collection during automatic washing procedure
- Ultra clear moulded blue overprinted graduations
- Anti suction base

Catalogue No	Cap., mL	O.D., mm	Height, mm	Subdivs., mL	Tolerance, mL	Pack qty
CYP-841-100B	100	39	160	2	2.0	5
CYP-841-250F	250	54	190	5	5.0	1
CYP-841-500Y	500	68	215	10	10.0	1
CYP-841-750Y	1,000	80	282	20	20.0	1

NOTE: Exposure to temperatures over 80°C may affect accuracy over time.



Fisherbrand®
Cylinders, polypropylene, ultra clear, spouted, hexagonal base, ISO6706

- With moulded graduations
- S.P.I. recycling code and date stamped for full traceability



Catalogue No	Cap., mL	Grads., mL	Pack qty
FB57950	10	0.2	10
FB57951	25	0.5	10
FB57952	50	1	10
FB57953	100	1	5
FB57954	250	2	each
FB57955	500	5	each
FB57956	1,000	10	each
FB57957	2,000	20	each



Fisherbrand®
Cylinders, PMP, transparent, with spout, hexagonal base

NEW

- With moulded graduations

Catalogue No	Cap., mL	Grads., mL
FB70350	10	0.2
FB70351	25	0.5
FB70352	50	1
FB70353	100	1
FB70354	250	2
FB70355	500	5
FB70356	1,000	10
FB70357	2,000	2