

# Flasks



Flasks, Pyrex® glass glass, flat bottom, narrow neck, BS2734, ISO1773

- With beaded rim

Catalogue No	Cap., mL	Neck I.D., mm	Height, mm	Pack qty
FHB-335-030L	100	18	110	5
FHB-335-050F	250	30	140	5
FHB-335-070W	500	30	170	2
FHB-335-090Q	1,000	36	200	2
FHB-335-110N	2,000	45	250	each
FHB-335-120K	6,000	57	340	each
FHB-335-130H	10,000	57	400	each



**Fisherbrand®**

Flasks, borosilicate glass, conical, narrow neck, ISO1773

- Erlenmeyer, graduated

Catalogue No	Cap., mL	Neck O.D., mm	Pack qty
FB33130	50	22	10
FB33131	100	22	10
FB33132	250	34	10
FB33133	500	34	10
FB33134	1,000	42	each
FB33135	2,000	50	each
FB33136	3,000	50	each
FB33137	5,000	50	each



Flasks, borosilicate glass, conical, wide neck

**Fisherbrand®**

- Erlenmeyer, graduated

Catalogue No	Cap., mL	Pack qty
FB33138	50	10
FB33139	100	10
FB33140	250	10
FB33141	500	10
FB33142	1,000	each
FB33143	2,000	each



Flasks, Pyrex® glass, conical, narrow neck, 1130 series, BS2734

- Erlenmeyer, graduated

Catalogue No	Cap., mL	Neck I.D., mm	Height, mm	Pack qty
FHB-375-070R	50	85	22	10
FHB-375-090L	100	105	22	10
FHB-375-130C	250	140	34	10
FHB-375-150T	500	175	34	10
FHB-375-170N	1,000	220	42	2



Flasks, Pyrex® glass, conical, wide neck, DIN 12385

- Erlenmeyer, graduated
- Supplied in packs of 10

Catalogue No	Cap., mL	Neck O.D., mm	Height, mm
FHB-485-030X	100	34	105
FHB-485-050Y	250	50	140
FHB-485-070C	500	50	175