



Gloves, vinyl, disposable, Suretouch™ PPE Category I

- Tear resistant beading for easy pull-on
- Special dispenser-style cases
- Ideal for food preparation, laboratory and pharmaceutical work
- Ambidextrous
- Supplied in packs of 100 gloves

| Catalogue No | Size |
|--------------------------|--------|
| Neutral | |
| SAR-760-010Y | Small |
| SAR-760-020F | Medium |
| SAR-760-030C | Large |
| Powder free inner coated | |
| SAR-765-010B | Small |
| SAR-765-020V | Medium |
| SAR-765-030S | Large |



Gloves, latex, disposable, Semperguard, PPE Category III

- Available in lightly powdered and powder free
- Rolled cuff, 0.34mm thickness at the palm and an AQL of 1.5
- Ambidextrous
- Chemical resistance data available on request
- Supplied in packs of 100 gloves

| Catalogue No | Size |
|--------------------------|-------------|
| Lightly powdered | |
| SAR-244-010P | Extra small |
| SAR-244-030J | Small |
| SAR-244-050D | Medium |
| SAR-244-070U | Large |
| SAR-244-090X | Extra Large |
| Powder free inner coated | |
| SAR-246-010D | Extra small |
| SAR-246-030U | Small |
| SAR-246-050X | Medium |
| SAR-246-070Y | Large |
| SAR-246-090C | Extra large |

Fax: 01509 231893

Telephone: 01509 233344

Gloves, nitrile, disposable, powder free, purple, PPE Category III




- Latex and protein free
- Textured fingertips for improved grip
- Available in 250mm length and 300mm length
- 0.10mm thickness
- Exceeds the requirements of EN374
- Three times better puncture resistance than latex
- Supplied in packs of 100 gloves



| Catalogue No | Size |
|--------------|-------------------|
| 250mm length | |
| FB69260 | Extra extra Small |
| FB69261 | Extra Small |
| FB69262 | Small |
| FB69263 | Medium |
| FB69264 | Large |
| FB69265 | Extra Large |
| 300mm length | |
| FB69290 | Extra extra Small |
| FB69291 | Extra Small |
| FB69292 | Small |
| FB69293 | Medium |
| FB69294 | Large |
| FB69295 | Extra Large |

Gloves, nitrile, disposable, non-sterile Touch 'n' Tuff™, PPE Category III

- Available in powdered and powder free
- Provide three times the puncture resistance of natural latex and vinyl gloves
- No additives to contaminate product or irritate skin
- Rolled cuff, 0.1mm thickness at the palm
- Supplied in packs of 100 gloves



| Catalogue No | Size |
|--|-------------|
| Powder free, green | |
| SAR-480-100R | Small |
| SAR-480-200N | Medium |
| SAR-480-300J | Large |
| SAR-480-400F | Extra large |
| Powdered, green | |
| SAR-480-010S | Small |
| SAR-480-020P | Medium |
| SAR-480-030M | Large |
| SAR-480-040J | Extra large |
| Powder free, blue with textured fingertips | |
| SAR-605-100E | Small |
| SAR-605-200A | Medium |
| SAR-605-300T | Large |
| SAR-605-400P | Extra large |