

Spatulas, metal

Fisherbrand®

- All spatulas have 3.6mm diameter except FB65069 which has 3mm diameter

Catalogue No	Length, mm	Width, mm	Pack qty
Nickel, drawn out flat to form narrow, bent blade at one end ❶			
FB65018	165	-	5
Stainless steel, lightweight, flat with rounded ends ❷			
FB65019	90	20	5
FB65021	150	30	5
Stainless steel, heavyweight, with one spatulate and one spoon end ❸			
FB65023	120	18	5
FB65026	150	18	5
FB65027	210	18	each
Nickel, Chattaway ❹			
FB65029	100	10	5
FB65031	150	10	5
FB65033	180	10	5
FB65035	200	10	5
FB65036	230	10	5
FB65098	270	10	5
Stainless steel, Chattaway ❺			
FB65037	100	10	5
FB65039	150	10	5
FB65041	180	10	5
FB65071	200	10	5
FB65043	230	10	5
FB65099	300	10	5
FB65100	350	10	5
FB65069	150	8	5
FB65061	100	4	5
FB65081	176	20	5
Stainless steel, with scoop and narrow bent blade			
FB65079	140	10	10
Stainless steel with offset spoon and chisel type blade			
FB65083	190	20	each



Spatulas, PTFE coated, stainless steel

Fisherbrand®

NEW

- Ideal for use under severe laboratory conditions
- Chemically resistant and non-contaminating
- Excellent non-wetting properties for easy cleaning and can be steam or chemically sterilised
- PTFE coated spatulas and tongs, providing a non-stick surface to which very few solid or liquid substances will permanently adhere



Catalogue No	Description
FB65110	150mm Chattaway spatula light
FB65111	200mm Chattaway spatula light
FB65112	100mm Chattaway spatula medium
FB65113	150mm Chattaway spatula medium
FB65114	180mm Chattaway spatula medium
FB65115	200mm Chattaway spatula medium
FB65116	230mm Chattaway spatula medium
FB65117	350mm Chattaway spatula medium
FB65118	100mm Chattaway spatula micro
FB65119	150mm Chattaway spatula micro
FB65120	200mm Chattaway spatula micro
FB65121	150mm trough spatula
FB65122	175mm trulla spatula
FB65123	120mm spoon one end
FB65124	150mm spoon one end
FB65125	210mm spoon one end
FB65126	100mm micro spoon
FB65127	150mm micro spoon
FB65128	150mm straight tongs
FB65129	200mm straight tongs
FB65130	150mm bow tongs
FB65131	200mm bow tongs