



Magnetic followers, glass

- Supplied in packs of 5

Catalogue No	Length, mm	Dia., mm
FB55624	6	5
FB55625	12	5
FB55626	19	5
FB55627	25	6
FB55628	40	8
FB55629	50	8



Magnetic followers, PTFE, centre rim

- Supplied in packs of 5, except FB55677 and FB55648 supplied individually

Catalogue No	Length, mm	Dia., mm
FB55646	50	8
FB55645	45	8
FB55648	70	10
FB55647	60	10
FB55644	40	8
FB55641	20	6
FB55640	12	4.5
FB55643	35	6
FB55642	25	6



Magnetic followers, PTFE, octagonal

- All have pivot rings, Alnico-V magnets encased in PTFE
- Supplied in packs of 10

Catalogue No	Length, mm	Dia., mm
FB55585	13	3
FB55587	13	8
FB55589	15	8
FB55591	25	8
FB55593	38	8
FB55595	51	8
FB55597	25	9.5
FB55599	38	9.5
FB55601	51	9.5
FB55603	64	9.5



Magnetic followers, PTFE, triangular

- Supplied in packs of 5

Catalogue No	Length, mm	Dia., mm
FB55650	12	6
FB55651	25	8
FB55652	35	9
FB55653	50	12



Magnetic followers, PTFE, single sided star head

- Supplied in packs of 5, except FB55664 and FB55665 supplied individually

Catalogue No	Length, mm	Dia., mm
FB55660	8	10
FB55661	12	14
FB55662	13	17
FB55663	15	25
FB55664	17	40
FB55665	17	60



Magnetic followers, PTFE, micro

- For stirring in test tubes, small size beakers or other containers
- Supplied in re-useable snap cap vials
- Supplied in packs of 10

Catalogue No	Length, mm	Dia., mm
FB55605	2	2
FB55606	5	2
FB55607	8	1.5
FB55609	7	2
FB55611	10	3



Magnetic followers, PTFE, value range

Catalogue No	Length, mm	Dia., mm	Pack qty
FB55930	12	3	10
FB55932	25	8	10
FB55934	40	8	10
FB55936	50	8	10