

# Precision Portable balances

## Scout® SKX Series



### The Next Generation of Portable Balance for your Classroom

- Superior Overload Protection & Practical Stackable Design for Convenient Storage
- Bright Backlit LCD & 4 Button Operation Make the New Scout Easy to See & Use
- Multiple connectivity options are available for data logging and communication
- Fast Weighing Speed and High Resolution Delivers Repeatable and Reliable Results

**Application** Weighing, Density Determination, Mole Weighing

**Display** Backlit LCD: 6-digit 7-segment with white LED backlight

**Operation** AC Adapter (Included), batteries (not included)

**Communication** Optional RS232, USB host, USB device, Bluetooth or Ethernet

**Construction** ABS housing, stainless steel pan, replaceable in-use cover

**Design Features** Fast stabilization time, menu lock switch, stackable design, security slot, integral weigh below hook for below balance weighing applications, transportation lock, removable stainless steel weighing platform, stability indicator, overload and underload indicators, low battery indicator, auto shut-off



Model	Capacity	Readability (d)	Platform	Test weight	Item №
SKX123	120g	0,001g	Ø 93 mm	100g, M1, included with balance	30253070
SKX222	220g	0,01g	Ø 120 mm	200g, M1, included with balance	30253071
SKX422	420g	0,01g	Ø 120 mm	200g, M1, included with balance	30253072
SKX622	620g	0,01g	Ø 120 mm	300g, M1, included with balance	30253073
SKX1202	1200g	0,01g	170mm x 140mm	1000g, M1	30253074
SKX2202	2200g	0,01g	170mm x 140mm	2000g, M1	30253075
SKX421	420g	0,1g	Ø 120 mm	200g, M1, included with balance	30253076
SKX621	620g	0,1g	170mm x 140mm	300g, M1, included with balance	30253077
SKX2201	2200g	0,1g	170mm x 140mm	2000g, M1	30253078
SKX6201	6200g	0,1g	170mm x 140mm	5000g, M1	30253079
SKX8200	8200g	1g	170mm x 140mm	8000g, M1	30253080

<b>Accessories</b>	<b>Item Nº</b>
RS232 Interface Kit, Scout	30268982
USB Host Interface Kit, Scout	30268983
USB Device Interface Kit, Scout	30268984
Bluetooth Interface Kit, Scout	30268985
Ethernet Interface Kit, Scout	30268986
Stacking Kit, set (6), Scout	30268987
Stacking Kit, set (1), Scout	30268988
Full-featured Impact Printer SF40A	30064202
SF40A Paper roll (57,5mm 2pcs)	12120799
SF40A Ink ribbon cassette	12120798
Auxiliary Display Kit, Scout	30269019
Security device (Laptop lock)	80850043
Top loading kit for Density Determination	30269020
Carrying case, Scout	30269021
In-use cover, Scout	30269022



SF40A IMPACT PRINTER