

Bags, biohazard



33



Heavy duty bags that can be loaded with used pipettes, Petri dishes, culture plates and membrane filters.

The entire load is then sterilised at a maximum of 121°C and incinerated. All bags are 32µm thick.

Red, large black biohazard symbol and precautionary procedures in 4 languages (English, French, German, Spanish) with sterilisation indicator patch which turns black at 121°C

Catalogue No	Width, mm	Height, mm	Pack qty
BAJ-530-010H	360	480	200
BAJ-530-020E	480	580	200
BAJ-530-030B	640	890	200

Clear, unprinted can be used as conventional bin liners

Catalogue No	Width, mm	Height, mm	Pack qty
BAJ-510-010V	300	610	100
BAJ-510-020S	610	760	100
BAJ-510-030P	610	910	100

Clear, standard biohazard warning symbol and precautionary procedures in 4 languages (English, French, German and Spanish)

Catalogue No	Width, mm	Height, mm	Pack qty
BAJ-520-010X	300	610	100
BAJ-520-020L	610	760	100
BAJ-520-030Y	610	910	100

Bag holders, biohazard, Clavies®



02



Recommended for use with the clear, unprinted and standard biohazard bag ranges (see above), these biohazard bag holders ensure the bag is securely supported on a stainless steel wire frame and polypropylene bottom. The bottom also serves as a containment tray in the event of a bag leakage. The tray is orange for easy identification as a biohazard containment system and can be autoclaved at 121°C (250°F) to maintain sterility. For ease of transportation they are supplied unassembled, however assembly only takes a few minutes and once done becomes permanent.

Catalogue No	Tray dimensions, mm	Height, mm	Bag size, mm
BAJ-510-500E	185 x 356	470	300 x 610
BAJ-510-505R	356 x 356	546	610 x 760
BAJ-510-510B	356 x 356	696	610 x 910

Bags, polyethylene with attached document pocket, minigrip

02



A range of specimen pouches featuring a secure, resealable section for the specimen and an open pocket for the associated documents.

Available either plain or printed with a biohazard warning. Opening along width. 50µm thickness.

Plain style

Catalogue No	Width, mm	Depth, mm	Pack qty
BAJ-350-090Y	150	140	100

Printed with biohazard warning

Catalogue No	Width, mm	Depth, mm	Pack qty
BAJ-350-070X	140	140	100
BAJ-350-110F	150	140	100

Bags, polyethylene, biohazard, B-safe

04



Ideal for heavy duty biohazard disposal.

Supplied in handy, bench top dispenser boxes.

Catalogue No	Width, mm	Depth, mm	Pack qty
BAJ-560-010J	255	355	50
BAJ-560-030D	455	760	50
BAJ-560-050U	790	1,015	25

Bags, polyethylene, biohazard, B-safe, Bin box

04



Ideal for heavy duty biohazard disposal.

Each bag is box shaped and has a double string pull. Supplied in handy, bench top dispenser boxes.

Catalogue No	Width, mm	Depth, mm	Height, mm	Pack qty
BAJ-560-070X	150	150	150	25
BAJ-560-090Y	200	200	250	25
BAJ-560-110F	490	270	400	25
BAJ-560-130W	280	280	510	25