



### Aluminium foil

**a** 03

Suitable for many laboratory applications, including sealing bottles, culture tubes, stopping leaks, reducing evaporation or heat losses, covering corks, making weighing boats or dishes.

Ideal for use in autoclaves.

Catalogue No	Length, m	Width, mm	Thickness, µm
<b>AKL-300-020P</b>	75	300	12
<b>AKL-300-030M</b>	75	450	12
<b>AKL-300-040J</b>	75	500	12
<b>AKL-300-045W</b>	75	600	12
<b>AKL-300-060D</b>	5	450	12
<b>AKL-310-040C</b>	75	100	18



### Aluminium foil, Super-Polyfoil™



**a** 03

For suspension cultures, media preparation or storage.

May be used to cover vessels for autoclaving and incubating.

Resists tearing and puncture, provides sterility protection, autoclavable and provides a chemically inert barrier. Polyester complies with FDA regulations governing direct food contact applications. A blue thermoplastic polyester layer is laminated to 1µm aluminium foil. The laminate is easier to remove following autoclaving than plain aluminium foil which is subject to oxidation.

Catalogue No	Width, mm	Length, mm	Pack qty
<b>AKL-400-010K</b>	102	102	100
<b>AKL-400-030E</b>	130	130	100
<b>AKL-400-050V</b>	203	203	100

### Flasks, shaker, polycarbonate



**s a NEW**

46

For suspension cultures, media preparation or storage

- Moulded in graduations
- Polycarbonate with polypropylene closure; vented closure has 0.22µm pore PTFE
- Sterile, non-pyrogenic and non-cytotoxic
- Certified sterile at 10<sup>6</sup> SAL USP Class VI (5 year shelf life)
- For single use; individually packed for easy storage and handling
- Available with plain bottom for use on bench top, or baffled bottom for shaker table use
- Available with vented or non-vented screw closures

Flask and non-vented closures are autoclavable.

Plain bottomed flask

Catalogue No	Capacity, mL	Style	Pack qty
<b>FB74440</b>	125	Vented	24
<b>FB74441</b>	250	Vented	12
<b>FB74442</b>	500	Vented	12
<b>FB74443</b>	1,000	Vented	6
<b>FB74444</b>	2,000	Vented	4
<b>FB74445</b>	2,800	Vented	4
<b>FB74446</b>	125	Non-vented	24
<b>FB74447</b>	250	Non-vented	12
<b>FB74448</b>	500	Non-vented	12
<b>FB74449</b>	1,000	Non-Vented	6

Baffled bottomed flask

Catalogue No	Capacity, mL	Style	Pack qty
<b>FB74450</b>	125	Vented	24
<b>FB74451</b>	250	Vented	12
<b>FB74452</b>	500	Vented	12
<b>FB74453</b>	1,000	Vented	6
<b>FB74454</b>	2,000	Vented	4