



### Disposable flasks, Erlenmeyer, PETG, sterile, Thermo Scientific Nalgene®

**Thermo**  
SCIENTIFIC

**NALGENE**



54

Optically clear PETG plastic. Ideal for shaker and suspension cell culture, media preparation and storage. Gamma sterilised, moulded-in graduations. Supplied with HDPE screw closures.

#### Plain bottom, white closure

| Catalogue No        | Alt. No   | Capacity, mL | Graduations, mL | Pack qty |
|---------------------|-----------|--------------|-----------------|----------|
| <b>TKV-137-005E</b> | 4112-0125 | 125          | 25              | 24       |
| <b>TKV-137-010L</b> | 4112-0250 | 250          | 50              | 12       |
| <b>TKV-137-020Y</b> | 4112-0500 | 500          | 50              | 12       |
| <b>TKV-137-025V</b> | 4112-1000 | 1,000        | 50              | 6        |
| <b>TKV-137-030F</b> | 4112-2000 | 2,000        | 500             | 4        |
| <b>TKV-137-030W</b> | 4112-2800 | 2,800        | 250             | 4        |

#### Baffled bottom, white closure

| Catalogue No        | Alt. No   | Capacity, mL | Graduations, mL | Pack qty |
|---------------------|-----------|--------------|-----------------|----------|
| <b>TKV-137-035S</b> | 4112-0125 | 125          | 25              | 24       |
| <b>TKV-137-040C</b> | 4113-0250 | 250          | 50              | 12       |
| <b>TKV-137-050W</b> | 4113-0500 | 500          | 50              | 12       |
| <b>TKV-137-055M</b> | 4113-1000 | 1,000        | 50              | 6        |
| <b>TKV-137-060T</b> | 4113-2000 | 2,000        | 500             | 4        |
| <b>TKV-137-065J</b> | 4113-2800 | 2,800        | 250             | 4        |

#### Plain bottom, blue vented closure

| Catalogue No        | Alt. No   | Capacity, mL | Graduations, mL | Pack qty |
|---------------------|-----------|--------------|-----------------|----------|
| <b>TKV-137-200G</b> | 4115-0125 | 125          | 25              | 24       |
| <b>TKV-137-210D</b> | 4115-0250 | 250          | 50              | 12       |
| <b>TKV-137-220A</b> | 4115-0500 | 500          | 50              | 12       |
| <b>TKV-137-230U</b> | 4115-1000 | 1,000        | 50              | 6        |
| <b>TKV-137-240R</b> | 4115-2000 | 2,000        | 500             | 4        |
| <b>TKV-137-245H</b> | 4115-2800 | 2,800        | 250             | 4        |

#### Baffled bottom, blue vented closure

| Catalogue No        | Alt. No   | Capacity, mL | Graduations, mL | Pack qty |
|---------------------|-----------|--------------|-----------------|----------|
| <b>TKV-137-300C</b> | 4116-0125 | 125          | 25              | 24       |
| <b>TKV-137-310W</b> | 4116-0250 | 250          | 50              | 12       |
| <b>TKV-137-320T</b> | 4116-0500 | 500          | 50              | 12       |
| <b>TKV-137-330Q</b> | 4116-1000 | 1,000        | 50              | 6        |
| <b>TKV-137-340N</b> | 4116-2000 | 2,000        | 200             | 4        |
| <b>TKV-137-350K</b> | 4116-2800 | 2,800        | 250             | 4        |

#### Vented closure, 0.2mm hydrophobic PTFE membrane

| Catalogue No        | Alt. No   | Description                        | Pack qty |
|---------------------|-----------|------------------------------------|----------|
| <b>TKV-137-100K</b> | 4114-0038 | Fits 250mL flasks                  | 12       |
| <b>TKV-137-110H</b> | 4114-0045 | Fits 500, 1,000 and 2,000mL flasks | 12       |



### Flasks, borosilicate glass, conical, culture, narrow neck

**Fisherbrand®**

41

Fluted, with four indents down the length of the wall.

| Catalogue No   | Capacity, mL |
|----------------|--------------|
| <b>FB55901</b> | 250          |
| <b>FB55903</b> | 500          |
| <b>FB55905</b> | 1,000        |
| <b>FB55907</b> | 2,000        |