



Bottles, borosilicate glass, laboratory, reagent and media, crimp top



50

Ideal for reaction vessels, sample storage and freeze drying applications.

Conforms to EP Type 1 standards.

Clear or amber bottles with crimp top and a wide range of seals and stoppers. Autoclavable.

Clear

Catalogue No	Alt. No	Capacity, mL	I.D. x O.D., mm	Dimensions [dia. x h], mm	Pack qty
BTF-780-050F	223738	5	13 x 20 (neck)	23 x 47	288
BTF-780-070W	223739	10	13 x 20 (neck)	25 x 54	288
BTF-780-090Q	223742	20	13 x 20 (neck)	32 x 58	288
BTF-780-110N	223743	30	13 x 20 (neck)	37 x 63	288
BTF-780-130H	223745	50	13 x 20 (neck)	43 x 73	288
BTF-780-150B	223747	100	13 x 20 (neck)	52 x 95	144
BTF-780-170S	223748	125	13 x 20 (neck)	54 x 107	144

Amber

Catalogue No	Alt. No	Capacity, mL	I.D. x O.D., mm	Dimensions [dia. x h], mm	Pack qty
BTF-782-050Q	223760	5	13 x 20 (neck)	23 x 47	288
BTF-782-070K	223761	10	13 x 20 (neck)	25 x 54	288
BTF-782-090E	223762	20	13 x 20 (neck)	32 x 58	288
BTF-782-110B	223763	30	13 x 20 (neck)	27 x 63	288
BTF-782-130S	223764	50	13 x 20 (neck)	43 x 73	288
BTF-782-150M	223766	100	13 x 20 (neck)	52 x 95	144

Rubber stoppers - straight plug style

Catalogue No	Alt. No	Description	For neck dimensions, mm	Pack qty
BTF-784-050E	224100-172	Natural red rubber	13 x 20 (I.D. x O.D.)	1,000
BTF-784-070V	224100-173	Grey chlorobutyl/45	13 x 20 (I.D. x O.D.)	1,000
BTF-784-090P	224100-175	Teflon® coated grey chlorobutyl-isoprene blend/40	13 x 20 (I.D. x O.D.)	1,000

Rubber stoppers - 2 leg lyophilisation

Catalogue No	Alt. No	Description	For neck dimensions, mm	Pack qty
BTF-784-110M	224100-192	Grey chlorobutyl-isoprene blend/50	13 x 20 (I.D. x O.D.)	1,000

Rubber stoppers - 3 leg lyophilisation

Catalogue No	Alt. No	Description	For neck dimensions, mm	Pack qty
BTF-784-130G	224100-202	Grey chlorobutyl/50	13 x 20 (I.D. x O.D.)	1,000

Aluminium crimp seals

Catalogue No	Alt. No	Description	For neck dimensions, mm	Pack qty
BTF-785-050V	24178-01	Open top, natural	20 (O.D.)	1,000
BTF-785-070P	224183-01	Centre disc tear out, natural	20 (O.D.)	1,000
BTF-785-090J	224191	Solid top, natural	20 (O.D.)	1,000
BTF-785-110G	224193-01	Tear off, natural	20 (O.D.)	1,000



For crimpers and decappers, please see page 739