

Cell/Tissue Culture and Cell Transfection

Pipettes - Serological pipettes, plastic



Pipettes, serological, plugged, polystyrene, paper peel



56

- Printed with product code, lot number and expiry date
- Paper backing allows for easy access and disposal after use
- Sterile by irradiation

Paper peel packaging, plugged, individually wrapped

Catalogue No	Alt. No	Description	Negative graduations, mL	Inner pack qty	Pack qty
PMK-335-011D	40501K	1mL pipette	-0.2	100	1,000
PMK-335-015S	40502K	2mL pipette	-0.3	100	1,000
PMK-335-030W	47505	5mL pipette	-2.0	50	200
PMK-420-040C	47510	10mL pipette	-3.0	50	500
PMK-335-050Q	47525	25mL pipette	-10.0	25	200
PMK-335-060N	47550	50mL pipette	-5.0	25	50



Pipettes, serological, disposable glass, Pyrex®

CORNING



02

- Calibrated to deliver with blow-out
- All pipettes have negative graduations
- Steriview wrappers (paper/plastic)
- Colours conform to ISO standards
- Long, slender, tapered tips to make pipetting go faster and easier

Individually wrapped, plugged, sterile

Catalogue No	Alt. No	Capacity, mL	Sleeve	Pack qty
TKV-600-030F	7077-5N	5	1	120
TKV-600-040C	7077-10N	10	1	100

Multi-plug, plugged, sterile

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Inner pack qty	Pack qty
TKV-610-080G	7078-1N	1	0.1	50	1,000
TKV-610-070J	7078-1CN	1	0.01	50	1,000
TKV-610-060M	7078-2N	2	0.01	35	700
TKV-610-090D	7078-5N	5	0.1	30	960
TKV-610-050P	7078-10N	10	0.1	20	120

Pipettes, serological, disposable, glass, Shorty

CORNING



66

Ideal for use in confined areas such as hoods.

Sterile, plugged pipettes with tooled tops.

Negative graduations make it applicable for both serological and measuring procedures. I.S.O. color coding is printed directly on Steri-View (one side paper, one side plastic) wrapper.

Individually wrapped

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Negative graduations, mL	Pack qty
TKV-605-010E	7077B-1	1	2/10	-4	200
TKV-605-020B	7077B-10	10	2/10	-4	200
TKV-605-045Y	7077B-25	25	2/10	-5	100
TKV-605-050P	7077B-50	50	5/10	-3	25