



### Pipette fillers, Thermo Scientific Matrix CellMate

**Thermo**  
SCIENTIFIC

50

- Separate aspirate and dispense controls
- Lithium ion battery technology
- Backlit LCD display
- 1mL pipette stabiliser
- Customisable ID area

Includes table stand, wall mount, 1mL pipette stabiliser, universal charging transformer.

#### Technical Specification - General

Aspirate speed, s	8
Dispensing speed, s	8 (+ gravity dispense)
Battery	Lithium ion
Run time, hr	20 (continuous use between charges approx.)
Recharge time, min	60 (approx)
Weight, g	220
Nosepiece Pipette Holder	Autoclavable silicone rubber
Pipette Compatibility	All major brands of glass or plastic serological pipettes, 1mL to 100mL
Filter assembly, µm	0.45
Fill time, s	5 (for 50mL pipette)

Catalogue No	Alt. No	Colour
<b>PMX-170-010Y</b>	9501	White
<b>PMX-170-020F</b>	9511	Clear
<b>PMX-170-030C</b>	9521	Blue
<b>PMX-170-040W</b>	9531	Red
<b>PMX-170-050T</b>	9541	Green

#### Accessories

Catalogue No	Alt. No	Description
<b>PMX-174-010H</b>	9057	Hydrophobic filters, 0.45µm, pack of 24
<b>PMX-174-020E</b>	9064	Nosepiece holder, silicone gripper, filter
<b>PMX-174-030B</b>	9065	Pipette gripper, silicone, pack of 4
<b>PMX-174-040V</b>	9066	Wall mount
<b>PMX-174-050S</b>	9067	Stand
<b>PMX-174-060P</b>	9068	Power supply universal, 100 to 240V 50/60Hz
<b>PMX-174-070M</b>	9070	Pipette support 1mL



### Pipette fillers, PipetBoy Acu

**INTEGRA**

**S**

04

Cable free, pipette controller.

Extremely lightweight, compact and easy to use. Works on a rechargeable 9V NiMH battery which is easily inserted at the top of the handle. One charge is sufficient to fill and purge at least 500 x 50mL pipettes. Accepts standard pipettes ranging from 0.1mL to 100mL. Comes in a range of five transparent colours and a traditional light grey to aid identification. With battery, transformer for 230V/50Hz, sterile filter, and wall mount.

Catalogue No	Alt. No	Colour
<b>PMR-100-010C</b>	155000	Light grey
<b>PMR-100-030T</b>	155015	Green
<b>PMR-100-050N</b>	155016	Red
<b>PMR-100-070H</b>	155017	Blue
<b>PMR-100-090B</b>	155018	Clear
<b>PMR-100-095X</b>	155019	Purple

#### Accessories

Catalogue No	Alt. No	Description
<b>PMR-100-510F</b>	153005	Rechargeable battery 9V NiMH
<b>PMR-100-520C</b>	151005	Handle mount
<b>PMR-100-530W</b>	151020	Pipettor, rubber, 0.1mL to 50mL
<b>PMR-100-115L</b>	153015	Filter, sterile, for PipetBoy

Please contact Fisher Scientific to discuss your EU compatible product requirements.