

Cell/Tissue Culture and Cell Transfection

Apparatus and consumables - TC-treated plates

Microplates, 96 well and 384 well, polystyrene, solid FluoroNunc™/LumiNunc™, Nunc™ Δ



S 76

- Plate and lid designed to reduce evaporation
- Minimum cross-talk
- White plates provide maximum reflection, minimum autofluorescence and autoluminescence
- Black plates provide minimum background in fluorescence reading and minimum back scattered light
- TC-treated



96 well, 128mm x 86mm, flat bottom, polystyrene, sterile

Catalogue No	Alt. No	Lid	Total well volume, mL	Volume, working, µL	Culture area, cm ²	Colour	Inner pack qty	Pack qty
TKT-186-010C	136101	Yes	400	50 to 250	0.33	White	1	50
TKT-186-020W	136102	Yes	400	50 to 250	0.33	White	10	160
TKT-186-030T	137101	Yes	400	50 to 250	0.33	Black	1	50
TKT-186-040Q	137103	Yes	400	50 to 250	0.33	Black	10	160

384 well, 128mm x 86mm, polystyrene, sterile

Catalogue No	Alt. No	Lid	Total well volume, mL	Volume, working, µL	Colour	Inner pack qty	Pack qty
TKT-183-010U	164610	Yes	120	10 to 100	White	10	30
TKT-183-030X	164564	Yes	120	10 to 100	Black	10	30
TKT-183-040L	165195	No	120	10 to 100	White	25	100
TKT-183-060F	164571	No	120	10 to 100	Black	25	100

Microplates, 96 and 384 well, polystyrene/polymer base, optical bottom, coated



76

- Ready to use Poly-D-Lysine or Collagen I pre-coated plate
- Consistent growth surface lot-to-lot
- High quality surface promotes cell attachment
- Produced in clean environment
- Optically clear plates with polymer bottom

Plates have a shelf life of seven months at 4°C to 25°C. Storage at 4°C to 8°C is recommended.

96 well

Catalogue No	Alt. No	Colour	Surface	Inner pack qty	Pack qty
PDL-152-036E	152039	Clear	Poly-D-Lysine	5	20
PDL-152-037C	152037	Black	Poly-D-Lysine	5	20
PDL-152-038A	152028	White	Poly-D-Lysine	5	20
COL-152-035S	152038	Clear	Collagen I	5	20
COL-152-036Q	152036	Black	Collagen I	5	20

384 well

Catalogue No	Alt. No	Colour	Surface	Inner pack qty	Pack qty
PDL-152-040N	152029	Black	Poly-D-Lysine	5	20
COL-152-040C	152041	Black	Collagen I	5	20



Microplates, 384 well, polystyrene, with optical bottom, Nunc™ Δ



S 76

- TC-treated
- Working volume 10 to 100µL/well
- Rounded square wells eliminate wicking

White and black upper structures with optically clear bottoms. Can be read from the bottom, top or both. Polystyrene 384 well upper structure and polymer base.

Catalogue No	Alt. No	Colour	Lid	Inner pack qty	Pack qty
TKT-184-010X	142762	White	Yes	10	30
TKT-184-020L	142761	Black	Yes	10	30
MPA-550-050E	164730	Black	No	10	30

