


Cell/Tissue Culture and Cell Transfection

Apparatus and consumables - TC-treated slides



Microscope slides, plastic, TC-treated, Thermo Scientific Nunc®



76

- Surface treatment ensures optimal conditions for cell attachment and growth
- Dimensions 25mm x 75mm

Catalogue No	Alt. No	Description	Inner pack qty	Pack qty
TKT-210-700Y	160005	Permanox™	20	100
TKT-210-710F	160004	Polystyrene	20	100



Chamber slides, TC-treated, Thermo Scientific Nunc®, Lab-Tek™, Chamber Slide™



82

Developed for microscopy of in-situ tissue culture.

- Cells grow on a microscope slide
- Growth chamber can be removed when culturing is complete
- Useful for viral and mycoplasma testing, chromosome studies, toxicity tests and immunocytology
- Certified surface treatment ensures optimal conditions for cell attachment and growth
- Permanox™ and glass do not interfere with fluorescent labels



TPX = Permanox™
PS = Polystyrene

Catalogue No	Alt. No	No. of wells	Slide material	Volume, working, mL	Culture area, cm²	Inner pack qty	Pack qty
TKT-210-430L	177372	1	Glass	2.0 to 4.5	9.4/well	8	96
CBR-177-305Y	177372SP	1	Glass	2.0 to 4.5	9.4/well	8	16
TKT-210-510N	177410	1	TPX	2.5 to 4.5	9.4/well	8	96
CBR-177-310P	177410SP	1	TPX	2.5 to 4.5	9.4/well	8	16
TKT-210-450F	177380	2	Glass	1.2 to 2.0	4.2/well	8	96
CBR-177-380R	177380SP	2	Glass	1.2 to 2.0	4.2/well	8	16
TKT-210-530H	177429	2	TPX	1.2 to 2.0	4.2/well	8	96
CBR-210-530Y	177429SP	2	TPX	1.2 to 2.0	4.2/well	8	16
TKT-210-470W	177399	4	Glass	0.5 to 0.9	1.8/well	8	96
CBR-177-399T	177399SP	4	Glass	0.5 to 0.9	1.8/well	8	16
TKT-210-550B	177437	4	TPX	0.5 to 0.9	1.8/well	8	96
CBR-177-315F	177437SP	4	TPX	0.5 to 0.9	1.8/well	8	16
TKT-210-490Q	177402	8	Glass	0.2 to 0.4	0.8/well	8	96
TKT-210-816M	177402SP	8	Glass	0.2 to 0.4	0.8/well	8	16
TKT-210-570S	177445	8	TPX	0.2 to 0.4	0.8/well	8	96
CBR-210-180J	177445SP	8	TPX	0.02 to 0.4	0.8/well	8	16
TKT-210-630D	178599	16	Glass	0.1 to 0.2	0.4/well	8	96
CBR-177-320M	178599SP	16	Glass	0.1 to 0.2	0.4/well	8	16
TKT-210-590M	177453	Flaskette	Glass	2.5 to 5.0	10/well	8	96
CBR-177-325C	177453SP	Flaskette	Glass	2.5 to 5.0	10/well	8	16
TKT-120-010L	170920	SlideFlask	PS	2.5 to 5.0	9.0/well	5	50

Accessories

Catalogue No	Alt. No	Description	Inner pack qty	Pack qty
TKT-200-210J	171862	Cover slip for TKT-120-010L	100	1,000
TKT-210-210C	171080	Cover slip for TKT-210-630D	-	96

All products are for research use only unless otherwise indicated.