

Cell scrapers

Fisherbrand®

S

03

- Perfect for gently harvesting cells from tissue culture flasks, when enzymes cannot be used
- Available in two sizes
- With polyethylene blades that swivel for even contact with growth surface
- Medical grade polystyrene handles ensure compatibility with a wide variety of culture flasks
- Pre-sterilised
- Individually wrapped

Catalogue No	Handle length, mm	Blade length, mm	Pack qty
FB55198	390	30	100
FB55199	250	18	100



Cell lifter

Fisherbrand®

S

41

- Maintains even contact with growth area in culture dishes
- Allows cell layers to be removed quickly and without damage
- Disposable lifters made entirely of polypropylene
- Pre-sterilised
- Individually wrapped

Catalogue No	Description	Handle length, mm	Blade length, mm	Pack qty
FB55160	Cell lifter	180	20	100



Cell scrapers and lifter

CORNING

S

66

- Unique blade design minimises cell damage
- Useful when harvesting cells
- Even contact with the growth surface
- Large scraper is designed for larger flasks
- Individually wrapped
- Sterile (SAL 10⁻⁶, certified non-pyrogenic)
- Cell lifter is especially useful with dishes and cell culture clusters

Catalogue No	Alt. No	Description	Blade length, mm	Handle length, mm	Pack qty
TKV-127-511T	3008	Cell lifter	19	180	100
TKV-127-513P	3010	Small scraper	18	250	100
TKV-127-515L	3011	Large scraper	30	390	100



For centrifuges, please see the Fisher Scientific Laboratory Catalogue, or visit our website.

