

Purification kit, miRNA, low-throughput, Nucleospin miRNA



86

NucleoSpin® miRNA is suitable for the parallel isolation of small RNA, large RNA and protein.

Superior selectivity - purify your RNA fractionated by size

Isolation of small RNA only (18 to 200 bases)

Isolation of small (18 to 200 bases) and large (>200 bases) RNA in two fractions

Isolation of total RNA (small and large in one fraction)

Additional isolation of total protein

Denatured protein ready-to-use for analysis in SDS-PAGE and Western blot.

Excellent RNA recovery - even without phenol/chloroform

Lysis with chaotropic salt

Spin column based procedure

Convenient handling - complete equipment provided

NucleoSpin filters - efficient sample homogenisation

NucleoSpin protein removal columns - highest purity of small RNA

rDNase - highly efficient on-column digestion of genomic DNA at RT

NucleoSpin, for parallel isolation of small and large RNA

Catalogue No	Alt. No	Pack qty
NZ74097110	740971.10	10
NZ74097150	740971.50	50
NZ7409711250	740971.250	250

NucleoSpin RNA II, for isolation of total RNA. Including rDNase and shredders

Catalogue No	Alt. No	Pack qty
NZ74095510	740955.10	10
NZ74095550	740955.50	50
NZ740955250	740955.250	250

NucleoSpin RNA/protein, for simultaneous isolation of total RNA and protein from unsplit samples (including rDNase and shredders)

Catalogue No	Alt. No	Pack qty
NZ74093310	740933.10	10
NZ74093350	740933.50	50
NZ740933250	740933.250	250

NucleoSpin TriPrep, for simultaneous isolation of total RNA, genomic DNA, and total protein from a wide variety of unsplit samples

Catalogue No	Alt. No	Pack qty
NZ74096610	740966.10	10
NZ74096650	740966.50	50
NZ740966250	740966.250	250

NucleoSpin RNA XS, for isolation of highly concentrated total RNA from extremely small amounts of starting material (elution volume down to 5µL)

Catalogue No	Alt. No	Pack qty
NZ74090210	740902.10	10
NZ74090250	740902.50	50
NZ740902250	740902.250	250

NucleoSpin FFPE RNA, for isolation of total RNA from formalin fixed, paraffin embedded samples. Including Paraffin Dissolver, Decrosslinking Buffer, and rDNase.

Catalogue No	Alt. No	Pack qty
NZ74096910	740969.10	10
NZ74096950	NZ740969.50	50
NZ740969250	NZ740969.250	250

Protein quantification assay

Catalogue No	Alt. No	Pack qty
NZ74096750	740967.50	50
NZ740967250	740967.250	250