

### Zeba™ micro desalt spin columns, Thermo Scientific Pierce



PN

Fast, high recovery desalting for tiny samples.

Thermo Scientific Zeba™ micro desalt spin columns are specially designed microcentrifuge columns containing 75µL of Zeba Resin that quickly and effectively process small volume (2µL to 12µL) protein samples with excellent recovery.

- Multiple sample processing in less than six minutes
- For use with small sample volumes of 2µL to 12µL
- Minimal sample dilution

Even low concentration protein samples (25µg/mL) are efficiently desalted without significant sample loss. The micro desalt spin columns provide for exceptional desalting performance, including >95% retention of salts and other small molecules of less than 1,000MW. They are optimised for use with proteins and other macromolecules >7,000MW.

Catalogue No	Pack qty
<b>PN89877</b>	25
<b>PN89878</b>	50



For a full range of centrifuges, tubes and accessories, please see the Fisher Scientific Laboratory Catalogue, or visit our website.



For dialysis products, please see the Cell Extraction and Lysis section, pages 35 to 49.

