



Bulk pack



Racked open

### Pipettor tips, universal fit

**Fisherbrand®**

50

#### SureOne Rack

- Flush fit inserts sit flush with rack, allowing easy multichannel pipetting
- Notched clasp designed to prop lid up for autoclave venting
- Sure stack lid rim prevents racks from sliding around
- Safe lock latch allows for one-handed opening and closing ability and won't pop open when dropped
- Sure snap design securely holds insert in place, preventing it from being lifted out while pipetting
- Flip-lid can be either a hinged or lift off lid; hinged format does not prevent use of multichannel pipettor

#### SureOne Insert

- Contact clip prevents unnecessary contact with the tips while clipping in inserts
- Lot numbered inserts for traceability include tip volume and Fisher Catalogue No
- Indented insert holes inhibit tip movement, stabilising them for multichannel work
- Rigid inserts are all 96 well multichannel format offering a robust non-pliable design providing secure tip seal

#### SureOne Reload Packaging

- Dust flap prevents debris from contaminating the product
- Canyon cut package allows for easy tip insert removal

A comprehensive range of non-filter and aerosol barrier pipettor tips with the environmentally friendly reload system equalling less waste.

Pipettor tips have a universal fit and are compatible with most popular brands of pipettor.

#### Standard Sterility Claims

Non sterile products - certified to be free of RNase/DNase and DNA

Sterile products - e-beam sterile products are certified to be free of RNase/DNase, Pyrogen, Bioburden, PCR\* inhibitors, Endotoxin free.

#### Filter tips

##### Racked sterile

Catalogue No	Description	Volume, $\mu$ L	Colour	Inner pack qty	Pack qty
<b>FB78098</b>	Micropoint tip, graduated at 2 $\mu$ L	0.1 to 10	Clear	96	960
<b>FB78100</b>	Micropoint tip, graduated at 2.5 $\mu$ L, extended length	0.1 to 10	Clear	96	960
<b>FB78102</b>	Micropoint tip Eppendorf style	0.1 to 10	Clear	96	960
<b>FB78104</b>	Micropoint tip Eppendorf style	2 to 20	Clear	96	960
<b>FB78106</b>	Bevelled tip	2 to 20	Clear	96	960
<b>FB78108</b>	Bevelled tip	10 to 100	Clear	96	960
<b>FB78110</b>	Bevelled tip, graduated at 10, 50, 100 $\mu$ L	20 to 200	Clear	96	960
<b>FB78112</b>	Micropoint tip, graduated at 100, 200, 500, 1,000 $\mu$ L	100 to 1,000	Clear	96	960

#### Standard tips

##### Micropoint tip, graduated at 2 $\mu$ L. 0.1 $\mu$ L to 10 $\mu$ L

Catalogue No	Pack type	Colour	Inner pack qty	Pack qty
<b>FB78000</b>	Bulk bag	Clear	-	1,000
<b>FB78002</b>	Racked	Clear	96	960
<b>FB78004</b>	Racked sterile	Clear	96	960
<b>FB78006</b>	Reload system pack	Clear	96	960

##### Extended length Micropoint tip, graduated at 2.5 $\mu$ L. 0.1 $\mu$ L to 10 $\mu$ L

Catalogue No	Pack type	Colour	Inner pack qty	Pack qty
<b>FB78008</b>	Bulk bag	Clear	-	1,000
<b>FB78010</b>	Racked	Clear	96	960
<b>FB78012</b>	Racked sterile	Clear	96	960

##### Micropoint tip Eppendorf style. 0.1 $\mu$ L to 20 $\mu$ L

Catalogue No	Pack type	Colour	Inner pack qty	Pack qty
<b>FB78014</b>	Bulk bag	Clear	-	1,000
<b>FB78016</b>	Racked	Clear	96	960
<b>FB78018</b>	Racked sterile	Clear	96	960

\*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche