



Reservoirs, reagent

Fisherbrand®

S

41

Disposable reagent reservoirs for use with multichannel pipettors.

- Can hold up to 100mL
- Clear PVC or white polystyrene
- Available in sterile or non-sterile packs

Catalogue No	Description	Pack qty
FB60600	50mL, PVC, clear, non-sterile	100
FB60601	50mL, PS, sterile, individually packed	80
FB60602	50mL, PS, natural, sterile, inner packs of 5	200
FB60603	100mL, PS, white, sterile, inner packs of 5	200



Reservoirs, reagent, disposable

CORNING

S

66

Sterile, modified polystyrene reservoirs for use with multichannel pipettors.

Catalogue No	Alt. No	Capacity, mL	Colour	Pack qty
PMP-331-010C	4870	50	Natural	200
PMP-331-020W	4871	50	Natural	100
PMP-331-030T	4872	100	White	200
PMP-331-040Q	4873	100	White	100



Reservoirs, reagent

**Thermo
SCIENTIFIC**

a

04

Catalogue No	Alt. No	Description	Pack qty
PMP-210-050S	9510027	Reagent reservoir, 'V' shaped, autoclavable	5
PMP-210-070M	9510037	8 channel reservoir, autoclavable for 4, 8 or 12 channel pipettors	10
PMP-210-080J	9510047	Single channel reservoir, autoclavable	10