

Microplates, 96 well, polystyrene, non-treated, Thermo Scientific Nunc® MicroWell™



76

Used in many fields including in antibiotic screens, serological tests, storing and screening of compounds and DNA libraries.

- Pinchbar design
- Choice of well formats in clear polystyrene: F, U or V
- Excellent optical quality
- The wells have raised rims and lids (if included) are supplied with rings which help prevent cross contamination
- Lids designed for easy orientation and stacking

Without lid

Catalogue No	Alt. No	Well shape	Volume, total µL/well	Sterile, Yes/No	Inner pack qty	Pack qty
DIS-984-070S	249662	Conical	300	Yes	1	50
DIS-984-090M	269620	Flat	400	No	5	60
DIS-984-110J	269787	Flat	400	Yes	1	50
DIS-984-060V	249570	Conical	300	No	10	180
DIS-984-250Q	249940	Conical	300	Yes	10	80
DIS-984-115W	260210	U96	300	Yes	10	180
DIS-984-030H	262162	U96	300	Yes	1	50
DIS-984-020K	268152	U96	300	No	10	180

With lid

Catalogue No	Alt. No	Well shape	Volume, total µL/well	Sterile, Yes/No	Inner pack qty	Pack qty
DIS-984-230W	277143	Conical	300	Yes	1	50
DIS-984-080P	249952	Conical	300	Yes	10	160
DIS-984-240T	249935	Conical	300	Yes	10	80



Microplates, 96 well, polystyrene, assay, clear

CORNING



66

General assay plates, ideal for solution based assays (i.e. haemagglutination or precipitation assays), serial dilutions, and general storage applications.

- Sterile plates are sterilised by gamma radiation and certified non-pyrogenic
- Flat, round or U-shaped well bottom

DPS-100-040C and **DPS-100-070Q** are processed to improve hydrophilicity, e.g. haemagglutination assay. Plates are TC treated where specified.

Catalogue No	Alt. No	Well shape	Lid	Treated, [Y/ N]	Sterile, Yes/No	Volume, total µL/well	Inner pack qty	Pack qty
DPS-100-010L	3367	Round, U	No	No	Yes	330	1	50
DPS-110-060M	3591	Flat	No	No	No	360	1	50
DPS-100-020Y	3370	Flat	Yes	No	Yes	360	20	100
DPS-100-030F	3795	Round, U	No	No	Yes	330	20	100
DPS-100-040C	3798 [†]	Round, U	No	No	No	330	25	100
DPS-100-050W	3896	V bottom	No	No	Yes	320	1	48
DPS-100-060T	3897	V bottom	No	No	No	320	25	100
DPS-100-070Q	3898	V bottom	No	No	No	320	25	100
DPS-100-100K	3628	Flat	Yes	Yes	Yes	360	20	100
DPS-100-110H	3788	Round, U	Yes	No	Yes	330	20	100

[†]Processed to improve hydrophilicity for haemagglutination and similar assays.

