

Microplates, 96 well, optical bottom, non-treated, Thermo Scientific Nunc®



S 76

White and black upper structures with optically clear bottoms. Can be read from top, bottom or both. Polystyrene upper structure is fitted with no. 1.5 German glass or polystyrene film. 96 well footprint for compatibility with automated handling equipment. Alphanumeric well identification. Optimised for scintillation counting. Dimensions 128mm x 86mm. Total volume 400µL and recommended working volume up to 200µL. Plates are non-sterile. Supplied without lids.

Catalogue No	Alt. No	Colour	Base	Inner pack qty	Pack qty
MPA-540-010A	265302	White	Polystyrene	10	30
MPA-540-020U	265301	Black	Polystyrene	10	30
MPA-540-030R	265300	Black	Coverglass	5	30

Accessory lids

Catalogue No	Alt. No	Colour	Sterile, Yes/No	Inner pack qty	Pack qty
MPA-540-503A	255983	Clear	No	20	60




Microplates, 96 well, without lid, flat bottom, non-sterile

CORNING

66

- Thin optically clear bottom allows direct viewing on microscope, readings down to 340nm
- Dense walls, white optimal for luminescence, black for fluorescence
- Maximum well volume 360µL

Catalogue No	Alt. No	Pigment	Binding property	Inner pack qty	Pack qty
DPS-132-010B	3631	Black	Non-treated	25	100
DPS-132-020V	3632	White	Non-treated	25	100
DPS-132-030S	3601	Black	High bind	25	100
DPS-172-050H	3651	Black/clear bottom	NBS™	25	100
DPS-172-030N	3604	White/clear bottom	NBS™	25	100
DPS-184-030R	3340	Black	CellBIND®	5	50



Microplates, 384 well, polymer, optical bottom, non-treated, Thermo Scientific Nunc®



76

- White and black upper structures with optically clear bottoms
- Can be read from the top, bottom or both
- Polystyrene upper structure and film
- Dimensions 128mm x 86mm
- Non-sterile

Catalogue No	Alt. No	Colour	Inner pack qty	Pack qty
MPA-550-010Q	242763	White	10	30
MPA-550-020N	242764	Black	10	30



Lids

Catalogue No	Alt. No	Description	Colour	Sterile, Yes/No	Inner pack qty	Pack qty
DIS-976-501M	448143	Long side with cut off corners	White	No	20	120
DIS-976-507A	264611	Long side with condensation ring and evaporation barrier	Clear	Yes	1	60
DIS-976-509T	264612	Long side with evaporation barrier	Clear	No	20	180
MPA-740-504F	250002	For 384/1536 well plates	Clear	Yes	20	60
MPA-740-506B	250003	For 384/1536 well plates	Clear	No	20	60