# **Histology**

## Cytology



Single



Double



Mask

### Microscope slides, Thermo Scientific Shandon Cytoslides®

## Thermo



- All slides are 75mm x 25mm, made from high quality optical glass
- Specimen area (12mm diameter) is outlined for rapid location during screening
- Available as coated and un-coated slides
- Coated slides eliminate the need for albumin coating procedures and the background staining necessary on albuminised, IHC or ISH processed slides
- Uncoated slides are available with the circle on the front or back

For use with the Thermo Scientific Shandon Cytospin® cytocentrifuge.

### Single Cytoslides®

Specially made for use with EZ Cytofunnels. A white, 12mm diameter circle on the slide defines the area where cells will be displayed to allow for efficient and rapid location. Uncoated slides available with circle on front or back.

Specially designed for use with EZ Double Cytofunnel® chambers. Each specimen is displayed onto the slide within a circle marked '1' or '2' which corresponds with the EZ Double Cytofunnel® chamber.

### Mask Cytoslides®

Specially produced for immunofluorescence. Black coating masks out UV light so that only the specimen within the circle can be viewed through a microscope.

Catalogue No	Alt. No	Description	Treated, [Y/ N]	Pack qty
SD5991056	5991056	Cytoslides®, single circle, circle on back	Υ	100
SD5991051	5991051	Cytoslides®, single circle, circle on back	N	100
SD5991059	5991059	Cytoslides®, single circle, circle on front	N	100
SD5991055	5991055	Cytoslides®, double circle	Υ	100
SD5991054	5991054	Cytoslides®, double circle	N	100
SD5991057	5991057	Cytoslides®, single circle, black mask	Υ	100
SD5991050	5991050	Cytoslides®, double circle, black mask	Υ	100

### Collection fluid, cytology, Thermo Scientific Shandon Cytospin®

# Thermo





- Ready-to-use alcohol/CARBOWAX® (PEG) combination for direct specimen collection and fixation
- Permits storage of specimens or unstained slides for weeks without adverse effect
- Packaged in several forms, including convenient pump bottle and pre-filled specimen cups

Catalogue No	Alt. No	Pack type	Quantity	Pack qty
SD6768001	6768001	With pump	3.8L	1
SD6768315	6768315	-	500mL	2
SD9990310	9990310	-	10L	2
SD9990323	9990323	Cup, filled to 60mL	120mL	125

### Red collection fluid, CytoRich® Red, Thermo Scientific Shandon

# Thermo



Lyses RBCs and solubilises proteins.

- For a higher cell yield and a cleaner background
- Ready to use alcohol based fluid
- Compatible with immunohistochemistry
- Pre-filled containers available

Catalogue No	Alt. No	Description	Quantity	Pack qty
SDB9990800	B9990800	With pump	3.8L	1
SDB9990801	B9990801	-	500mL	2
SDB9990802	B9990802	Cups, filled to 10mL	20mL	180
SDB9990803	B9990803	Cups, filled to 60mL	120mL	125

All products are for research use only, unless otherwise