

### Microscope slide folder

**Fisherbrand®**

46

Cardboard slide trays for 76mm x 26mm slides. The horizontal compartments are recessed protecting the specimens when the cover is closed. 20 place.

Catalogue No	Dimensions, mm
<b>FB70901</b>	208 x 341 x 10



### Microscope slide boxes, plastic, colour coded

**Fisherbrand®**

46

Range of durable plastic microscope slide boxes ideal for storage or transport. Stackable.

Catalogue No	Description	Lining	Capacity, slides	Colour
<b>FB70860</b>	Push-fit lid	Cork	25	Blue
<b>FB70861</b>	Push-fit lid	Cork	25	Green
<b>FB70862</b>	Push-fit lid	Cork	25	Red
<b>FB70863</b>	Push-fit lid	Cork	25	White
<b>FB70866</b>	Thumb latch lock	Cork	50	Blue
<b>FB70867</b>	Thumb latch lock	Cork	50	Red
<b>FB70868</b>	Thumb latch lock	Cork	50	White
<b>FB70870</b>	Nickel plated clasp	Cork	100	Blue
<b>FB70871</b>	Nickel plated clasp	Cork	100	Green
<b>FB70872</b>	Nickel plated clasp	Cork	100	Red
<b>FB70873</b>	Nickel plated clasp	Cork	100	Yellow
<b>FB70874</b>	Nickel plated clasp	Cork	100	White
<b>FB70875</b>	Nickel plated clasp	Cork	100	Grey
<b>FB70876</b>	Nickel plated clasp	Cork	100	Black
<b>FB70880</b>	Hinged lid	Unlined	25	Blue
<b>FB70881</b>	Hinged lid	Unlined	25	Red
<b>FB70882</b>	Hinged lid	Unlined	25	White
<b>FB70883</b>	Hinged lid	Unlined	100	Blue
<b>FB70884</b>	Hinged lid	Unlined	100	Red
<b>FB70885</b>	Hinged lid	Unlined	100	White
<b>FB70895</b>	Nickel-plated clasp	Foam lined	100	Blue
<b>FB70896</b>	Nickel-plated clasp	Foam lined	100	Green
<b>FB70897</b>	Nickel-plated clasp	Foam lined	100	Red
<b>FB70898</b>	Nickel-plated clasp	Foam lined	100	Yellow
<b>FB70899</b>	Nickel-plated clasp	Foam lined	100	White

